



Reagent set, aluminum, 50 mL sample, ECR method

Artikel-Nr.: 2443400
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Für die Bestimmung von Aluminium mit der ECR-Methode. Für 25 ml Probe.

For determination of Aluminum by the Eriochrome Cyanine R (ECR) Method. Hach Method 8326. Sample Size: 50 mL Range: 0.002 - 0.220 mg/L. Reagent set contains ECR Reagent Powder Pillows, Hexamethylenetetramine Buffer Reagent Powder Pillows and ECR Masking Reagent Solution for approximately 100 tests.

Technische Daten

Anzahl Bestimmungen:	100
Messbereich:	0 - 0,220 mg/L Al
Messgerät:	DR/3000, DR/2000
Methode:	Eriochrome Cyanine R
Methode (Nummer):	8326