



Digital titrator cartridge, sodium hydroxide, 3.636 N

Artikel-Nr.:

1438001

CHF Preis:

Kontakt

Kein Versanddatum angezeigt

Gefahrstoffe



For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Anzahl Bestimmungen:	40 - 100
Messbereich:	100 - 1.000 mg/L CO ₂
Methode:	Phenolphthalein
Methode (Nummer):	8205
Parameter:	Carbon Dioxide
Plattform:	Digital Titrator
Titrant:	3.636 N NaOH