



EDTA digital titrator cartridge, 0.0716 M

Artikel-Nr.:

2081701

Gefahrstoffe

CHF Preis:

Kontakt

Kein Versanddatum angezeigt



EDTA standard solution. For the determination of iron by Digital Titrator titration. 0.0716 M. Digits x 0.1 = mg/L Fe in 50 mL

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Anzahl Bestimmungen:	40 - 100
Messbereich:	10 - 100 mg/L Fe
Methode:	EDTA
Parameter:	Iron
Plattform:	Digital Titrator
Titrant:	0.0716 M EDTA