



Reagent set, iron (total), LR, digital titrator

Artikel-Nr.: 2449200

Gefahrstoffe

VERALTETER ARTIKEL
Dieser Artikel ist nicht mehr verfügbar.



For determination of iron by EDTA Digital Titrator titration. Hach Method 8214. Range: 10 - 100 mg/L as Fe. Reagent Set Contains: Citrate Buffer Powder Pillows, Sodium Periodate Powder Pillows, Sulfosalicylic Acid Powder Pillows and 0.0716 M EDTA Digital Titrator Cartridge. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Anzahl Bestimmungen:	100
Messbereich:	10 - 100 mg/L Fe
Methode:	TitraVer Titration
Methode (Nummer):	8214
Plattform:	Digital Titrator
Spezielles Ausstattungsmerkmal:	Iron (Total) Reagent set, Digital Titrator, LR