



Be Right™



Titrationst Kartusche (digital), Silbernitrat, 1,128 N, 1 Stück

Artikel-Nr.: 1439701
CHF Preis: Kontakt
Verfügbar

Gefahrstoffe



Gebrauchsfertige Titrierlösung für die Chlorid-Bestimmung (Silbernitrat-Methode), 10 - 10000 mg/L Cl⁻.

For chloride determination by Digital Titrator argentometric titration. Hach method 8207. Digits = mg/L Cl⁻ in 50 mL sample. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Anzahl Bestimmungen:	40 - 100
Messbereich:	100 - 10000 mg/L Cl ⁻
Methode:	Silver Nitrate
Methode (Nummer):	8207
Parameter:	Chloride
Plattform:	Digital Titrator
Titrant:	1.128 N AgNO ₃