



ICE-PIC .5NTU CVD

Artikel-Nr.:	5222500
CHF Preis:	Kontakt
Kein Versanddatum angezeigt	

Can be used with the 1720C and D instruments. Verify calibration in less than 60 seconds. No loss of sample flow to complete a verification and no sample cell to recondition after verification.

Saves time

Verification is the fastest available, taking less than one minute.

Ease of use

Compact and lightweight, ready for use at any time.

Accuracy

±10% or better, and each device is supplied with a Certificate of Accuracy.

Repeatability

Exceptional repeatability performance to be tracked closely.

Affordability

One-time investment, no consumables. One-year warranty.

Technische Daten

Inhalt:	ICE-PIC, Anweisungen.
---------	-----------------------

Inhalt

ICE-PIC, Anweisungen.