



Package Piercer & Sampler

Artikel-Nr.: 29981

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Headspace oxygen analyzer for pharmaceutical vials or ampoules. Semi-automatic motorized device reduces operator variability. Custom-made vial and ampoule holders available. To be used with 510.

Accurate oxygen analysis in less than one minute

Proven system

Headspace volumes as small as 1 mL

Technische Daten

Genauigkeit:	$\pm 1\%$ of reading or $\pm 0.01\%$ volume O ₂ whichever is greater
Gewährleistung:	24 Monate
Inhalt:	Package sampler only.
	Flow chamber 32017, piercing needle assembly 280582 (ampoules) and /or 280583 (vials) must be ordered separately.
Messbereich:	0.05% - 21% volume O ₂
Spannungsversorgung (Volt):	110 - 240 V AC
Wiederholbarkeit:	r ₉₅ = $\pm 0.08\%$ O ₂

Inhalt

Package sampler only. Flow chamber 32017, piercing needle assembly 280582 (ampoules) and /or 280583 (vials) must be ordered separately.