



pH test kit, 5.5 - 8.5 pH, model 17F

Artikel-Nr.: 147006
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technische Daten

Anzahl Bestimmungen:	200
Kleinste Schrittweite:	0.1 pH
Koffertyp:	B
Messbereich:	5,6 - 8,4 pH
Methode:	Color disc/ Bromthymol blue
Modell:	17F
Parameter:	pH, mid range
Plattform:	Color Disc
Test-Kit Modell:	17F
Versandgewicht:	1