



Test kit, phosphate, total ortho-/meta-, model PO-24

Artikel-Nr.: 225001
 CHF Preis: Kontakt
 Kein Versanddatum angezeigt

Gefahrstoffe



Kit provides measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Kit provides rapid measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Technische Daten

Anzahl Bestimmungen:	100 klare
	100 trübe
Kleinste Schrittweite:	0.02
	0.1
	1
Koffertyp:	F
Messbereich:	Triple Range
Messbereich (2):	0 - 4 mg/L PO ₄
Messbereich (3):	0 - 40 mg/L PO ₄
Methode:	Color disc/Ascorbic acid
Modell:	PO-24
Parameter:	Phosphate
Plattform:	Color Disc
Standard-Messbereich:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Test-Kit Modell:	PO-24
Versandgewicht:	4,5