



Iron color disc test kit, model IR-18B

Artikel-Nr.:

146401

CHF Preis:

Kontakt

Verfügbar

Gefahrstoffe



Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.25
Koffertyp:	B
Messbereich:	High Single Range
Methode:	Color disc/ 1,10 Phenanthroline
Modell:	IR-18B
Parameter:	Iron
Plattform:	Color Disc
Standard-Messbereich:	0.25 - 7 mg/L Fe
Test-Kit Modell:	IR-18B
Typ:	Single Metric