



Test Kit, chlorine (free & total) & pH, model CN-67

Artikel-Nr.: 1411100
CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

One of Hach's most popular test kits

Includes color discs and a comparator for monitoring chlorine and pH

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anwendung:	Pool and Spa
Anzahl Bestimmungen:	300
Messbereich:	6,6 - 8,4 pH-Einheiten
Messbereich (2):	0 - 3,4 mg/L
Modell:	CN-67
Parameter:	Chlor (frei & gesamt) und pH
Plattform:	Color Disc