



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Artikel-Nr.: 223002
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Testet auf Chlor (0 - 3,4 mg/L), Härte (1 - 20 gpg), Eisen (0 - 4 mg/L) und pH-Wert (4 - 10).

Technische Daten

Anwendung:	Water Quality
Anzahl Bestimmungen:	100
Messbereich:	0 - 3,4 mg/L
Modell:	CN-39WR
Parameter:	Chlor, Härte, Eisen und pH
Plattform:	Color Disc