



Testkit, glycol visual, model EG-1

Artikel-Nr.:

2186400

CHF Preis:

Kontakt

Kein Versanddatum angezeigt

Gefahrstoffe



Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	25
Koffertyp:	G
Messbereich:	Vorhandensein oder Fehlen von Glykol, 0 - über 150 mg/L
Methode:	Visual/ Appearance of color
Parameter:	Glycol - Ethylene glycol in oil or water
Test-Kit Modell:	EG-1
Versandgewicht:	3,5