



Test kit cyanide, model CYN-3

Artikel-Nr.: 201002
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures. This kit provides a simple, low-cost method using a stepped color comparator.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.01
Koffertyp:	R
Messbereich:	0 - 0,3 mg/L CN ⁻
Methode:	Color disc/ Pyridine-Pyrazolone
Modell:	CYN-3
Parameter:	Cyanide
Plattform:	Color Disc
Standard-Messbereich:	0 - 0.3 mg/L CN ⁻
Test-Kit Modell:	CYN-3
Versandgewicht:	2,5