



Test kit, free chlorine, model CN-70F

Artikel-Nr.: 1454201
CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Quickly and accurately measure free chlorine levels in the field or in the lab

Technische Daten

Anzahl Bestimmungen:	200
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kleinste Schrittweite:	0.02
	0.1
Koffertyp:	D
Messbereich:	Dual Range
Messbereich (2):	0 - 3.4 mg/L Cl ₂
Methode:	Color disc/DPD
Parameter:	Chlor Frei
Plattform:	Color Disc
Test-Kit Modell:	CN-70F ²