



Sulfite test kit, model SU-5

Artikel-Nr.: 148002
CHF Preis: Kontakt
Verfügbar

Gefahrstoffe



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	1
	10
Koffertyp:	D
Messbereich:	1.0 - 20 mg/L & 10 - 200 mg/L SO_3^{2-}
Methode:	Drop count, titration Iodometric
Parameter:	Sulfite - As SO_3^{2-}
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	SU-5
Versandgewicht:	2