



Test kit, sulfite, model SU-DT

Artikel-Nr.: 2063300
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste mögliche Schritte:	0.4-4.0**
Koffertyp:	F
Messbereich:	4 - 800 mg/L SO_3^{2-}
Methode:	Digital Titrator/Iodometric
Parameter:	Sulfite - As SO_3^{2-}
Plattform:	Digital Titrator
Test-Kit Modell:	SU-DT