



Be Right™



Test kit, sulfate, model SF-1

Artikel-Nr.:

225100

Gefahrstoffe

CHF Preis:

Kontakt

Lieferung innerhalb von 2 Wochen



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	50
Koffertyp:	D
Messbereich:	50 - 200 mg/L
Methode:	Extinction/Turbidimetric
Parameter:	Sulfate - As SO ₄ ²⁻
Test-Kit Modell:	SF-1
Versandgewicht:	1,5