



ELECTRODE PKG, SODIUM ION

Artikel-Nr.:

2348100

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Sodium Analysis Package includes Sodium Platinum Series Combination Electrode, two ammonium chloride reference cartridges, 100 and 1000 mg/L standards, reagents, beaker, and manual.

Everything you need for ISE measurements in one package.

Hach's ISE Analysis Packages provide a complete solution for ion selective measurements, including ISE, standard solutions, ionic strength adjustor, and replacement electrolyte cartridges.

Maximum life through reduced maintenance.

Sension Platinum Series probes' replaceable gel reference electrolyte cartridge is a fast and easy maintenance solution that maximizes the life of your electrode.

Accurate

You don't just need an answer, but the right one, and fast.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Technische Daten

Diaphragma:	Open
Durchmesser:	12 mm
Elektroden-Typ:	Laboratory combination
Elektrolytlösung:	2595802
Gewährleistung:	6 Monate
ISA erforderlich:	4451569
Kabellänge:	1 m
Länge:	240 mm
Material:	Sensor Body: PMP
Messbereich:	0.01 - saturation mg/L Na ⁺
Sensortyp:	Glass
Spezielles Ausstattungsmerkmal:	Platinum Series
Temperaturbereich:	0 - 50 °C
Thermistor:	-

Erforderliches Zubehör

- Ammoniumchlorid-Referenzkartusche, 2 Stück (Item 2595802)
- Natrium/Kalium-Ionenstärke-Regler, Pulverkissen, 100 Stück (Item 4451569)