



PROBE, DISSOLVED OXYGEN, 1M CABL

Artikel-Nr.:

5197000

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Clark-Type Polarographic Dissolved Oxygen probe with 1 meter cable and 5-pin circular connector compatible with Hach Sension meters.

Accurate

You don't just need an answer, but the right one, and fast.

Optimized stabilization algorithms in Hach meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and the quickest response, even in challenging environments.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand in the future.

Technische Daten

Elektroden-Typ:	Laboratory Polarographic Dissolved Oxygen
Elektrolytlösung:	2759126
Gewährleistung:	6 Monate
Kabellänge:	1 m
Länge:	149 mm
Material:	Sensor Body: Epoxy
Messbereich:	0 - 20 mg/L
Sensortyp:	Platinum cathode, silver anode with replaceable PTFE membrane module (item no. 5197300)
Temperaturbereich:	0 - 50 °C

Erforderliches Zubehör

- Membranes (Item 5197300)
- Fülllösung, innen, polarographische Sonden für gelösten Sauerstoff, 50 mL (Item 2759126)