



Conductivity probe (3-m)

Artikel-Nr.:

5197503

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Graphite, 4-pole conductivity probe with 3 meter cable and 5-pin circular connector compatible with Hach Sension meters.

Accurate

You don't just need an answer, but the right one, and fast.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

Technische Daten

Durchmesser:	15 mm
Elektroden-Typ:	Laboratory
Gewährleistung:	6 Monate
Kabellänge:	3 m
Länge:	120 mm
Material:	Sensor Body: Noryl
Messbereich:	Conductivity: 0.1 µS/cm - 201 mS/cm
Messbereich (2):	TDS: 0 - 50000 mg/L as NaCl
Sensortyp:	4-pole graphite
Temperaturbereich:	N/A N/A N/A