



## Rucksack für tragbare Messgeräte, groß, Nylon, Räder, mit Fächern

Artikel-Nr.: 2990800

CHF Preis: Kontakt  
Lieferung innerhalb von 2 Wochen

Für den sicheren Transport von tragbaren Messgeräten, Zubehör und Standard-Lösungen.

### Bequem

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

### Individuell anpassbar

The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it.

### Ergonomisch

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

---

## Technische Daten

|                          |                          |
|--------------------------|--------------------------|
| Abmessungen (H x B x T): | 432 mm x 356 mm x 152 mm |
| Gewicht:                 | Capacity: 13,6 kg        |
| Länge:                   | 152 mm                   |
| Material:                | Nylon, plastic textile   |