



Reagent set, hardness (total), digital titrator, full range

Artikel-Nr.: 2272000

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Für die Bestimmung von Gesamt-Härte, 10-160 mg/l CaCO₃. EDTA Methode.

For Total Hardness determination by EDTA Digital Titrator titration with Calmagite indicator. Range: 10-160 and 100-4000 mg/L as CaCO₃. Set contains ManVer® 2 Indicator, Hardness 1 Buffer Solution and one cartridge each of 0.0800 M and 0.800 M EDTA.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Messbereich:	10 - 160 mg/L & 100-4000 mg/L CaCO ₃
Methode:	EDTA Titration
Parameter:	Härte (Total)
Spezielles Ausstattungsmerkmal:	Reagent set, Digital Titrator