



Analysenwaage DENVER SI-234, 230 g Wägebereich, 0,1 mg Ablesbarkeit

Artikel-Nr.:

SMSI-234

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Schalengröße: 90 mm, Datenausgang: RS 232

Linearity: 0.2 mg, Reproducibility: 0.1 mg, Calibration: internal, simply by pressing a key, Equipment: large, backlit display, easy to use large rectangular breeze shield, below weighing, bubble level, afj. feet, anti-theft security log, power supply. Reporting according to GMP/GLP, Applications: counting, formulation, target/percent-weighing, animal weighing, weighing in different weighing units.