



Spare parts kit 2 years, 8810 ISE Spare parts kit for 2 years, 8810 ISE NH4 analyser

Artikel-Nr.: 368810,93006

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Spare parts kit_2 years, 8810 ISE analyser w/ ammonia electrode composed of »1 x 359110,01650 » Stirrer complete 2 x 359090,80000 »
Tubing holder cassette 1 x 359110,01400 » Siphon tube for reactor vessel

Technische Daten

Zusammensetzung:

- 1 x 359110,01650 » Stirrer complete
- 2 x 359090,80000 » Tubing holder cassette
- 1 x 359110,01400 » Siphon tube for reactor vessel
- 1 x 359110,01500 » Overflow tube for reactor vessel
- 2 x 429=110=040 » VITON O-ring 11x 4 mm (around above)
- 2 x 359012,00070 » PTFE washer 18 x 12.5 x 3 mm (around above)
- 3 x 359015,20035 » Threaded sleeve PG 13.5 red
- 1 x 08810=A=0010 » Conn. 90°, G $\frac{1}{4}$ ", stem 6 mm tube, PP
- 1 x 587=008=002 » Conn. straight, G $\frac{1}{4}$ ", 6/8 mm tube, PP
- 1 x 359103,10070 » Conn. 90°, G $\frac{1}{8}$ ", 4/6 mm tube, PP
- 1 x 359103,10077 » Conn. 90°, G $\frac{3}{8}$ ", 10/12 mm tube, PP
- 1 x 587=114=003 » Conn. 90°, G $\frac{3}{8}$ ", stem 12 mm tube, PP
- 1 x 359103,10132 » Conn. straight, 1/8" MNPT, 6/8 mm tube, PP
- 1 x 359110,02300 » Conn. straight, G $\frac{3}{8}$ ", Stem 12 mm tube, PP
- 10 x 359500,11106 » Union fitting 2 x nipple 1.6 mm tube
- 4 x 359090,70015 » Tube set TYGON grey-grey (12 pieces)
- 3 x 151399,90002 » TYGON tubing 1,6 x 3,2 mm (per meter)
- 3 x 359110,04200 » Reagent injection tube, 8810 vessel

3 x 359110,04110 » FEP tubing 1,65 x 3 mm (per meter)

1 x 125=000=003 » Ammonia (NH₃) ISE electrode

1 x 125=000=103 » Internal electrolyte for ammonia electrode (1 liter)