



Be Right™



Nitrat-Stickstoff-Standardlösung, 1.000 mg/L NO₃-N (NIST), 500 mL

Artikel-Nr.:

1279249

Gefahrstoffe

CHF Preis:
Verfügbar

Kontakt



Ion Selective Electrode (ISE) standard solutions for laboratory analyses

Zur Überprüfung der Genauigkeit von Nitrat-Bestimmungen. 1.000 mg/L NO₃-N; 4.430 mg/L NO₃-.

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Technische Daten

Bemerkungen:	Observed shelf life after opening, following recommended storage and use: 3 months
Beschreibung:	Nitrat-Stickstoff-Standardlösung (NIST)
Gebindeeinstellungen:	Bottle
Konzentration:	1.000 mg/L NO ₃ -N
Lagerbedingungen:	10 - 25 °C (protect from light)
Parameter:	Nitrate
Plattform:	Lösung
Volumen / Packungsgröße:	500 mL