



Copper Standard Solution, 10 mg/L as Cu (NIST), 100 mL MDB

Artikel-Nr.: 12932
CHF Preis: Kontakt
Verfügbar

APA6000 Copper Standard Solutions, 100 mL, 10 mg/L (NIST)

Technische Daten

Beschreibung:	Copper Standard Solution as Cu (NIST)
Konzentration:	10 mg/L
Lagerbedingungen:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Messgerät:	APA6000
Parameter:	Hardness
Volumen:	100 mL MDB
Volumen / Packungsgröße:	100 mL