



Cyanuric acid test kit, model CY-3

Artikel-Nr.: 185102

CHF Preis: Kontakt

Kein Versanddatum angezeigt

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anwendung:	Pool and Spa
Anzahl Bestimmungen:	100
Kleinste Schrittweite:	20
Koffertyp:	D
Messbereich:	20 - 100 mg/L
Methode:	Turbidimetry
Modell:	CY-3
Parameter:	Cyanuric Acid
Test-Kit Modell:	CY-3