



## Test kit tannin/lignin

Artikel-Nr.: 193701  
CHF Preis: Kontakt  
Kein Versanddatum angezeigt

Gefahrstoffe



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

**Tannin and lignin are analyzed together and reported as mg/L tannic acid.**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.5
	5
Koffertyp:	R
Messbereich:	0 - 15 mg/L as Tannic Acid
Messbereich (2):	0 - 150 mg/L
Methode:	Color disc/Tyrosine
Modell:	Tannin/Lignin
Parameter:	Tannin/Lignin
Plattform:	Color Disc
Standard-Messbereich:	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L
Test-Kit Modell:	TA-3