



Test kit, acidity, model AC-DT

Artikel-Nr.: 2064000
 CHF Preis: Kontakt
 Kein Versanddatum angezeigt

Gefahrstoffe



Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Technische Daten

Anzahl Bestimmungen:	100
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Kleinste mögliche Schritte:	0.1-1.10
Koffertyp:	F
Messbereich:	10 - 4,000 mg/L CaCO ₃
Methode:	Digital Titrator/Sodium hydroxide
Parameter:	Acidity - (P or Total) as CaCO ₃
Plattform:	Digital Titrator
Test-Kit Modell:	AC-DT