



Test kit, chloride, model CD-51

Artikel-Nr.: 208601
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kleinste Schrittweite:	500
	5000
Koffertyp:	B
Messbereich:	500 - 10000 mg/L & 5000 - 100000 mg/L Cl ⁻
Methode:	Drop count titration/ Silver nitrate
Parameter:	Chloride - As Cl ⁻
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	CD-51
Versandgewicht:	1