



Cylinder, graduated, polymethylpentene, 250 mL $\pm$ 1.4 mL, 2.0 mL divisions

Artikel-Nr.: 217246
CHF Preis: Kontakt
Kein Versanddatum angezeigt

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated "To Contain/To Deliver" at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Technische Daten

Gehäusematerial:	Polymethylpentene
Kapazität:	250 mL
Messmethode:	To Deliver/To Contain
Typ:	Single Metric
Volumen / Packungsgröße:	1 each
Zertifizierung:	Class B