



Test kit, formaldehyde, model FM-1

Artikel-Nr.: 2183100
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.05%
	0.5%
Koffertyp:	D
Messbereich:	0,05 - 1% & 0,5 - 10%
Methode:	Drop count titration/ Thymolphthalein
Parameter:	Formaldehyde - As CH ₂ O
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	FM-1
Versandgewicht:	2,5