



Test kit, chromium, model CH-14

Artikel-Nr.: 222702
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the drop count titration method to determine chromium (high range), this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	5
	50
Koffertyp:	D
Messbereich:	5 - 100 mg/L & 50 - 1000 mg/L Cr ⁶⁺
Methode:	Drop count titration/ Thiosulfate
Parameter:	Chromium, high range - As hexavalent Cr
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	CH-14
Versandgewicht:	2