



Test kit, nitrogen , model NI-8

Artikel-Nr.: 224100
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.1
Koffertyp:	B
Messbereich:	0 - 3.0 mg/L NH ₃ -N
Methode:	Color disc/Nessler reagent
Parameter:	Nitrogen, Ammonia, mid range - For freshwater
Plattform:	Color Disc
Test-Kit Modell:	NI-8