



Reagent set, carbon dioxide, digital titrator

Artikel-Nr.: 2272700
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



For determination of Carbon Dioxide by Digital Titrator titration. Hach Method 8205. Range: 10 - 1000 mg/L as CO₂ Set includes: Phenolphthalein Pillows and two Sodium Hydroxide Digital Titrator Cartridges (0.3636 N and 3.3636 N). Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Technische Daten

Anzahl Bestimmungen:	40 - 100
Messbereich:	10 - 1.000 mg/L CO ₂
Methode (Nummer):	8205
Plattform:	Digital Titrator