



Agar R2A, dehydrated 500 g

Artikel-Nr.: 2281034
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Low nutrient medium for enumeration and cultivation of bacteria from potable water. Used for pour or spread plate and MF methods; described in the 16th edition of Standard Methods. Recommended for predicting nitrification more accurately in potable systems using chloramines as a secondary disinfectant.

Technische Daten

Media Type:	Agar R2A,
Region:	US
Volumen / Packungsgröße:	500 g