



Reagent powder pillows, silver 1, pk/50

Artikel-Nr.: 2293566

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Reagent for determination of Silver by the Cadion 2B Colorimetric Method. Hach Method 8120. Range: 0.02 - 0.70 mg/L Ag. Sample Size: 50 mL. Pack of 50 Pillows. Also available in reagent set 2296600.

Technische Daten

Anzahl Bestimmungen:	50
Methode:	Colorimetric
Plattform:	Pulverkissen