



Test kit, ascorbic acid, model ASC-1

Artikel-Nr.: 2308100

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Gefahrstoffe



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Fußnote:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Kleinste Schrittweite:	10
Koffertyp:	D
Messbereich:	10 - 200 mg/L
Methode:	Drop count titration/ Iodometric
Parameter:	Ascorbic Acid
Plattform:	Tropfenzähl-Titration
Test-Kit Modell:	ASC-1
Versandgewicht:	2,5