



Test kit, professional boiler treatment/boiler feed, model BT-DT

Artikel-Nr.: 2350600

VERALTETER ARTIKEL
Dieser Artikel ist nicht mehr verfügbar.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. This kit provides testing capability for alkalinity, chloride, hardness, sulfite, phosphate, and pH, with a wide choice of ranges.

Offers testing of treatment chemicals and feedwater parameters, with the accuracy demanded by water treatment professionals

The included Hach Digital Titrator simplifies titration procedures, saves time, and provides greater accuracy than drop count testing

Kit also includes color disc for simple, fast analysis

Technische Daten

Anwendung:	Boiler Water
Plattform:	Digital Titrator