



Filter, membrane, pore size 0.2 µm, diameter 25 mm, 100/pk

Artikel-Nr.: 2387200

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Filters are treated with an irgalan black process that renders the filters dark gray, eliminates natural autofluorescence and aids visibility. These polycarbonate screen membrane filters exhibit extremely low non-specific adsorption. The fluorescing dye used does not discolor the filter. Manufacturing provides exceptionally clean membranes, as free as possible of viable and nonviable microorganisms. Diameter: 25 mm. Pore size: 0.2 µm. Poretics.

Technische Daten

Durchmesser:	25 mm
Filter Typ:	Filter
Material:	Polycarbonate
Porengröße:	0,2 µm
Volumen / Packungsgröße:	100 /Packung