



Standard solution, sulfite, 15 mg/L as SO₃ (equivalent), 500 mL

Artikel-Nr.: 2408449

CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Sulfite equivalent standard solution. Prepared with stabilized Sodium Thiosulfate. Equivalent to 15 mg/L as SO₃. 500 mL Bottle

Technische Daten

Beschreibung:	Sulfite Equivalent Solution as SO ₄ ²⁻ (NIST)
Konzentration:	15 mg/L
Plattform:	Lösung