



Reagent set, molybdenum, LR, 50 mL

Artikel-Nr.: 2440200
CHF Preis:
Kein Versanddatum angezeigt

Kontakt

Gefahrstoffe



For determination of Low Range Molybdate Molybdenum by the Ternary Complex method (Hach method 8169). Range: 0.02 - 3.00 mg/L as Mo. Sample Size: 50 mL. Set contains: Molybdenum 1 Reagent Powder Pillows and Molybdenum 2 Reagent Solution for approximately 100 tests.

Technische Daten

Anzahl Bestimmungen:	100
Messbereich:	0,02 - 3,00 mg/L Mo
Messgerät:	DR/3000, DR/2000
Methode:	Ternary Complex
Methode (Nummer):	8169
Probenvolumen:	50 mL