



Test kit, phenols, model PL-1

Artikel-Nr.: 2483600
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



Utilizing the simple color disc method with a range of 0 - 4 mg/L, this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.2 mg/L
Koffertyp:	D
Messbereich:	0 - 4.0 mg/L
Methode:	Color disc/ 4-aminoantipyrine
Modell:	PL-1
Parameter:	Phenols
Plattform:	Color Disc
Standard-Messbereich:	0 - 4.0 mg/L
Test-Kit Modell:	PL-1