



Test Kit, AccuVac® Color Disc, Ozone, 0.05-1.50 mg/L

Artikel-Nr.: 2518050
CHF Preis: Kontakt
Kein Versanddatum angezeigt

Gefahrstoffe



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Technische Daten

Anzahl Bestimmungen:	25
Kleinste Schrittweite:	0.05
Koffertyp:	G
Messbereich:	0 - 1,50 mg/L O ₃
Methode:	Color disc/ Indigo trisulfonate
Modell:	NA
Parameter:	Ozone
Plattform:	Accuvac
Standard-Messbereich:	0 - 1.50 mg/L O ₃
Test-Kit Modell:	AccuVac® Kit
Versandgewicht:	4,1