



Bottle rack, plastic-coated, 15 compartments

Artikel-Nr.: 2551800

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Technische Daten

Spezielles Ausstattungsmerkmal:	Rack, 15 Bottle Compartments
Volumen / Packungsgröße:	1