



Standard solution, manganese, 10 mg/L as Mn (NIST), pk/20, 2 mL PourRite® Ampules

Artikel-Nr.:

2605820

CHF Preis:

Kontakt

Verfügbar

Gefahrstoffe



For accuracy checks in manganese determinations

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Technische Daten

Beschreibung:	Manganese Standard Solution as Mn (NIST)
Konzentration:	10 mg/L
Plattform:	Lösung
Volumen / Packungsgröße:	mL/20