



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Artikel-Nr.: 2618800

CHF Preis: Kontakt

Kein Versanddatum angezeigt

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Technische Daten

Durchmesser:	47 mm
Filter Typ:	Filter
Material:	Polyethersulfone
Packing:	200 St.
Porengröße:	0,2 µm
Volumen / Packungsgröße:	200 /pk