



Economy kit m-ColiBlue24[®], glass ampules, 200 test

Artikel-Nr.: 2779200

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

In 24 hours, Hach's m-ColiBlue24[®] Broth gives you simultaneous results for total coliforms and E. coli. Use m-ColiBlue24[®] Broth to accurately monitor and evaluate total coliform and E. coli in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

With ready-to-use m-ColiBlue24[®] Broth, you can differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms. Minimize background growth of non-coliform bacteria.

Read and confirm results in 24 hours

Enumerate total coliforms and E. coli on one petri dish

Conduct emergency testing and routine monitoring

Get superior sensitivity - 1 CFU/100mL

Eliminate special equipment and ultraviolet lamp