



APA6000 Alkalinity Standard 2, 500 mg/L, 1L

Artikel-Nr.: 2826253

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

The APA6000 Alkalinity Process analyzer uses five reagents each of which is supplied in a 1L bottle designed to provide reagents for 30 days of continuous use. APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample. APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration. APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration. The fifth reagent, APA Alkalinity Cleaning Solution (2697053), cleans the system during the Prime and Instrument Clean functions.

Technische Daten

Konzentration:	500 mg/L
Lagerbedingungen:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Messgerät:	APA6000, EZ4004
Parameter:	Alkalinity